
Product Datasheet

**S6 Ribosomal Protein (phospho-Ser240/244)
antibody (Biotin) (orb501621)**

Description

S6 Ribosomal Protein (phospho-

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ELISA, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human RPS6 around the phosphorylation site of Ser240/244
Target	Phospho-RPS6 (Ser240 + Ser244)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	29kDa
Uniprot ID	P62753
Dilution Range	WB=1:100-1000 ,IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:100-1000
Expiration Date	12 months from date of receipt.