


Product Datasheet

Rat IL-6/Interleukin-6 ELISA Kit (orb50165)

| | |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | Rat IL-6/Interleukin-6 ELISA Kit |
| Reactivity | Rat |
| Range | 62.5 pg/ml - 4,000 pg/ml |
| Tested Applications | ELISA |
| Immunogen | Expression system for standard: E.coli; Immunogen sequence: F25-T211 |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| Note | For research use only |
| Application notes | Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 62.5 pg/ml |
| Sample Types | cell culture supernatants, serum and plasma (heparin, EDTA, citrate). |
| MW | 24357 MW |
| Uniprot ID | P20607 |
| Sensitivity | 5 pg/ml |
| Expiration Date | Please enquire. |



Rat IL-6
PicoKine
ELISA Kit
standard
cur...