

Product Datasheet

phospho-FOXO4 Thr451 antibody (Biotin) (orb501690)

Description	phospho-FOXO4 Thr451 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Biotin
Tested Applications	ICC, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human FOXO4 around the phosphorylation site of Thr451
Target	phospho-FOXO4 (Thr451)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P98177
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.