

Product Datasheet

PDCD4 (phospho-Ser457) antibody (Biotin) (orb501810)

Description	PDCD4 (phospho-Ser457) antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Biotin
Tested Applications	ICC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human PDCD4 around the phosphorylation site of Ser457
Target	phospho-PDCD4 (Ser457)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q53EL6
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.