

Product Datasheet

Vimentin (phospho-Ser55) antibody (Biotin) (orb501816)

Description	Vimentin (phospho-Ser55) antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Vimentin around the phosphorylation site of Ser55
Target	phospho-Vimentin (Ser55)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
CAS Number	Vimentin
Uniprot ID	P08670
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.