

---

## Product Datasheet

**LIM Kinase 2 (phospho-Thr505) antibody (Biotin)  
(orb501898)**

**Description**

Rabbit polyclonal antibody against

conjugated to Biotin

**Species/Host**

Rabbit

**Reactivity**

Rat

**Conjugation**

Biotin

**Tested Applications**

IHC-P

**Immunogen**

KLH conjugated synthesised phosphopeptide derived from human LIM Kinase 1 around the phosphorylation site of Thr505

**Target**

Limk2

**Preservatives**

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

**Form/Appearance**

Lyophilized or Liquid

**Concentration**

- 100 µg (in 100 µl): 1 mg/ml

**Storage**

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**

70kDa

**Uniprot ID**
**P53671**
**Dilution Range**

WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, ICC: 1:50-200, IF: 1:50-200

**Expiration Date**

12 months from date of receipt.