


Product Datasheet

Human ErbB-2 ELISA Kit (orb50195)

Description	Human ErbB-2 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: T23-T652
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:10, concentration in serum and plasma is 3107-9304 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA)
MW	137910 MW
Uniprot ID	P04626
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human ErbB-2 PicoKine ELISA Kit standard...