


Product Datasheet

Human Fetuin A/AHSG ELISA Kit (orb50196)

Description	Human Fetuin A/AHSG ELISA Kit
Reactivity	Human
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A19-V367
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:100000, concentration in serum and plasma is 219-805 ug/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	39325 MW
Uniprot ID	P02765
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human Fetuin A PicoKine ELISA Kit standa...