


Product Datasheet

Rat MIP-3 Alpha/CCL20 ELISA Kit (orb50207)

Description	Rat MIP-3 Alpha/CCL20 ELISA Kit
Reactivity	Rat
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S27-M96
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	10875 MW
Uniprot ID	P97884
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Rat MIP-3
alpha/CCL20
PicoKine
ELISA Kit...