


Product Datasheet

Mouse Rantes ELISA Kit (orb50208)

Description	Mouse Rantes ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S24-S91
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is 56-192 pg/ml
Sample Types	cell culture supernatants. serum and plasma (heparin, EDTA, citrate).
MW	10071 MW
Uniprot ID	P30882
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Mouse
Rantes
PicoKine
ELISA Kit
standard...