


Product Datasheet

Mouse RBP4/Retinol Binding Protein 4 ELISA Kit (orb50217)

Description	Mouse RBP4/Retinol Binding Protein 4 ELISA Kit
Reactivity	Mouse
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E19-L201
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:500, concentration in serum and plasma is 1-2 ug/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
MW	28405 MW
Uniprot ID	H7BWY6
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse
RBP4
PicoKine
ELISA Kit
standard
C...