

Product Datasheet

NPM (phospho-Thr199) antibody (Biotin) (orb502207)

Description	Rabbit polyclonal antibody against NPM (phospho-Thr199) conjugated to Biotin
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Biotin
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human NPM around the phosphorylation site of Thr199 RD(p-T)PA
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	38kDa
Uniprot ID	P06748
Dilution Range	WB: 1:50-200, ELISA: 1:100-1000
Expiration Date	12 months from date of receipt.