
Product Datasheet

Tuberin (phospho-Ser939) antibody (Biotin) (orb502252)

Description

Rabbit polyclonal antibody against

to Biotin

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Biotin

Tested Applications

ELISA, IHC-Fr, IHC-P, WB

Immunogen

KLH conjugated Synthesised phosphopeptide derived from human Tuberin around the phosphorylation site of Ser939

Target

TSC2

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

200kDa

Uniprot ID
P49815
Dilution Range

WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000

Expiration Date

12 months from date of receipt.