

Product Datasheet

MEF2C (phospho-Thr20) antibody (Biotin) (orb502415)

Description	MEF2C (phospho-Thr20) antibody (Biotin)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Biotin
Tested Applications	ELISA, IHC-Fr, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human MEF2C around the phosphorylation site of Thr20
Target	phospho-MEF2C (Thr20)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q06413
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.