

Product Datasheet

Caveolin 2 (phospho-Tyr19) antibody (Biotin) (orb502493)

Description	Caveolin 2 (phospho-Tyr19) antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from mouse Caveolin2 around the phosphorylation site of Tyr19.
Target	phospho-Caveolin-2 (Tyr19)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P51636
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.