

Product Datasheet

phospho-GluR1 Ser863 antibody (Biotin) (orb502540)

Description	phospho-GluR1 Ser863 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ELISA, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human GluR1 around the phosphorylation site of Ser863
Target	phospho-GluR1 (Ser863)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P42261
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.