

Product Datasheet

ERK1/2 antibody (Biotin) (orb502555)

| | |
|----------------------------|--|
| Description | ERK1/2 antibody (Biotin) |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Biotin |
| Tested Applications | ELISA, IHC-Fr, IHC-P |
| Immunogen | KLH conjugated Synthesised phosphopeptide derived from rat MAPK3 around the phosphorylation site of Tyr204 |
| Target | phospho-MAPK3 (Tyr204) |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q63844 |
| Dilution Range | IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000 |
| Expiration Date | 12 months from date of receipt. |