

Product Datasheet

PPAR Gamma (Phospho-ser273) antibody (FITC) (orb502604)

Description	PPAR Gamma (Phospho-ser273) antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human PPAR Gamma around the phosphorylation site of ser273
Target	Phospho-PPAR Gamma (ser273)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P37231
Dilution Range	IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.