

Product Datasheet

phospho-C-REL Ser503 antibody (FITC) (orb502639)

Description	Rabbit polyclonal antibody against phospho-C-REL Ser503 conjugated to FITC
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Tested Applications	ICC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human C-REL around the phosphorylation site of Ser503 TS(p-S)DS
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	69kDa
Uniprot ID	Q04864
Dilution Range	ICC: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.