

## Product Datasheet

# IRS1 (phospho-Ser323) antibody (FITC) (orb502688)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | IRS1 (phospho-Ser323) antibody (FITC)   |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human, Mouse  |
| <b>Conjugation</b>         | FITC  |
| <b>Tested Applications</b> | FC, IF  |
| <b>Immunogen</b>           | KLH conjugated synthesised phosphopeptide derived from human IRS1 around the phosphorylation site of Ser323 |
| <b>Target</b>              | phospho-IRS1 (Ser323)   |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.                                       |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Concentration</b>       | 1mg/ml  |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.                            |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Uniprot ID</b>          | <b>P35568</b>   |
| <b>Dilution Range</b>      | Flow-Cyt=1µg/Test,IF=1:100-500  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |