

---

## Product Datasheet

**PDCD4 (phospho-Ser67) antibody (FITC)  
(orb502706)**

**Description**

Rabbit polyclonal antibody against

FITC

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse

**Conjugation**

FITC

**Tested Applications**

ICC

**Immunogen**

KLH conjugated synthesised phosphopeptide derived from human PDCD4 around the phosphorylation site of Ser67

**Target**

PDCD4

**Preservatives**

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

**Form/Appearance**

Lyophilized or Liquid

**Concentration**

- 100 µg (in 100 µl): 1 mg/ml

**Storage**

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**

51kDa

**Uniprot ID**
[Q53EL6](#)
**Dilution Range**

ICC=1:50-200 IF=1:50-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.

**Expiration Date**

12 months from date of receipt.