

## Product Datasheet

# ERK1/2 (phospho-Thr183/185) antibody (FITC) (orb502887)

<b>Description</b>	ERK1/2 (phospho-Thr183/185) antibody (FITC)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	ICC, IF
<b>Immunogen</b>	KLH conjugated Synthesised phosphopeptide derived from mouse ERK1 around the phosphorylation site of Thr183/Tyr185
<b>Target</b>	phospho-ERK1 + 2 (Thr183/Tyr185)
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">Q63844</a>
<b>Dilution Range</b>	IF=1:100-500,ICC=1:100
<b>Expiration Date</b>	12 months from date of receipt.