

---

## Product Datasheet

**phospho-ATRIP Ser68+Ser72 antibody (FITC)  
(orb502901)**

**Description**

Rabbit polyclonal antibody against

conjugated to FITC

**Species/Host**

Rabbit

**Conjugation**

FITC

**Tested Applications**

ICC, IF

**Immunogen**

KLH conjugated synthesised phosphopeptide derived from human ATRIP around the phosphorylation site of Ser68 + Ser72 [LA(p-S)QAL(p-S)QC]

**Target**

ATRIP

**Preservatives**

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

**Form/Appearance**

Lyophilized or Liquid

**Concentration**

- 100 µg (in 100 µl): 1 mg/ml

**Storage**

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**

86kDa

**Uniprot ID**
**Q8WXE1**
**Dilution Range**

IF=1:100-500, ICC=1:100-500

**Expiration Date**

12 months from date of receipt.