
Product Datasheet

**Gria2 (phospho-Ser880) antibody (HRP)
(orb503130)**

Description

Rabbit polyclonal antibody against

HRP

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

HRP

Tested Applications

IHC-P, WB

Immunogen

KLH conjugated Synthesised phosphopeptide derived from human Gria2 around the phosphorylation site of Ser880

Target

GRIA2

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

96kDa

Uniprot ID
P42262
Dilution Range

WB=

Expiration Date

12 months from date of receipt.