

---

## Product Datasheet

**PLC gamma 2 (phospho-Tyr1217) antibody (HRP)  
(orb503333)**

**Description**

Rabbit polyclonal antibody against

conjugated to HRP

**Species/Host**

Rabbit

**Reactivity**

Human

**Conjugation**

HRP

**Tested Applications**

ELISA, IHC-Fr, IHC-P, WB

**Immunogen**

KLH conjugated Synthesised phosphopeptide derived from human PLC gamma 2 around the phosphorylation site of Tyr1217

**Target**

PLCG2

**Preservatives**

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

**Form/Appearance**

Lyophilized or Liquid

**Concentration**

- 100 µg (in 100 µl): 1 mg/ml

**Storage**

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**

139kDa

**Uniprot ID****P16885****Dilution Range**

WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000

**Expiration Date**

12 months from date of receipt.