

Product Datasheet

NFKB1 (phospho-Ser927) antibody (HRP) (orb503355)

Description	NFKB1 (phospho-Ser927) antibody (HRP)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	HRP
Tested Applications	ELISA, IHC-Fr, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from rat NFKB1 around the phosphorylation site of Ser927
Target	Phospho-NFKB1 (Ser927)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19838
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.