
Product Datasheet

**MAPK14 (phospho-Thr180/Tyr182) antibody (HRP)
(orb503533)**

Description

MAPK14 (phospho-Thr180/Tyr182)

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

HRP

Tested Applications

IHC-Fr, IHC-P, WB

Immunogen

KLH conjugated Synthesised phosphopeptide derived from human p38 MAPK around the phosphorylation site of Thr180/Tyr182

Target

Phospho-P38 MAPK (Thr180 + Tyr182)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID[Q16539](#)**Dilution Range**

WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500

Expiration Date

12 months from date of receipt.