

Product Datasheet

IKB alpha (phospho-Ser36) antibody (HRP) (orb503676)

Description	IKB alpha (phospho-Ser36) antibody (HRP)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	HRP
Tested Applications	ICC, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human IKB alpha around the phosphorylation site of Ser36
Target	phospho-IKB alpha (Ser36)
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P25963
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.