

---

## Product Datasheet

### Survivin (phospho-Thr34) antibody (HRP) (orb503707)

**Description**

Rabbit polyclonal antibody against

to HRP

**Species/Host**

Rabbit

**Conjugation**

HRP

**Tested Applications**

ELISA

**Immunogen**

KLH conjugated synthesised phosphopeptide derived from human Survivin around the phosphorylation site of Thr34

**Target**

BIRC5

**Preservatives**

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

**Form/Appearance**

Lyophilized or Liquid

**Concentration**

- 100 µg (in 100 µl): 1 mg/ml

**Storage**

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**

16kDa

**Uniprot ID**
[O15392](#)
**Expiration Date**

12 months from date of receipt.