

Product Datasheet

Smad3 (phospho-Ser425) antibody (HRP) (orb503725)

Description	Smad3 (phospho-Ser425) antibody (HRP)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	HRP
Tested Applications	ICC, IHC-Fr, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Smad3 around the phosphorylation site of Ser425
Target	Phospho-Smad3 (Ser425)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P84022
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.