
Product Datasheet

**CD18 (Phospho-Thr758) antibody (HRP)
(orb503855)**

Description

Rabbit polyclonal antibody against

	HRP
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human CD18 around the phosphorylation site of Thr758 SA(p-T)TT
Target	ITGB2
Preservatives	Lyophilized or Liquid
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	82kDa
Uniprot ID	P05107
Dilution Range	ELISA=1:5000-10000, IHC-F=1:100-500, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.