

Product Datasheet

ATF2 (Phospho-Thr69 + Thr51) antibody (HRP) (orb503960)

Description	ATF2 (Phospho-Thr69 + Thr51) antibody (HRP)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	HRP
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human ATF2 around the phosphorylation site of Thr69 + Thr51
Target	phospho-ATF2 (Thr69 / Thr51)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P15336
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.