
Product Datasheet

**p53BP1 (phospho-Ser25/29) antibody (HRP)
(orb504031)**

Description

p53BP1 (phospho-Ser25/29) antibody

| | |
|----------------------------|---|
| Species/Host | Rabbit |
| Reactivity | Mouse, Rat |
| Conjugation | HRP |
| Tested Applications | ELISA, ICC, IHC-Fr, IHC-P |
| Immunogen | KLH conjugated Synthesised phosphopeptide derived from human p53BP1 around the phosphorylation site of Ser25/29 |
| Target | phospho-p53BP1 (Ser25 + Ser29) |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 213kDa |
| Uniprot ID | Q12888 |
| Dilution Range | IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500,ELISA=1:5000-10000 |
| Expiration Date | 12 months from date of receipt. |