
Product Datasheet

pPP2A alpha antibody (PE) (orb504536)

Description

Rabbit polyclonal antibody against pPP2A

Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	PE
Tested Applications	IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human PPP2CA around the phosphorylation site of Tyr307
Target	PPP2CA
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	34kDa
Uniprot ID	P62714
Expiration Date	12 months from date of receipt.