

## Product Datasheet

# NFkB p65 (phospho-Thr435) antibody (PE) (orb504578)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | NFkB p65 (phospho-Thr435) antibody (PE)  |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human, Mouse   |
| <b>Conjugation</b>         | PE   |
| <b>Tested Applications</b> | FC, ICC, IF  |
| <b>Immunogen</b>           | KLH conjugated Synthesised phosphopeptide derived from human NFkBp65 around the phosphorylation site of Thr435 |
| <b>Target</b>              | phospho-NFkB p65 (Thr435)  |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.  |
| <b>Form/Appearance</b>     | Liquid   |
| <b>Concentration</b>       | 1mg/ml   |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.                               |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Uniprot ID</b>          | <a href="#">Q04206</a>   |
| <b>Dilution Range</b>      | Flow-Cyt=2ug/Test,IF=1:100-500,ICC=1:100   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |