

Product Datasheet

HNF4 (phospho-Ser313) antibody (PE) (orb504794)

Description	HNF4 (phospho-Ser313) antibody (PE)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE
Tested Applications	FC, ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human HNF4 alpha around the phosphorylation site of Ser313
Target	phospho-HNF4 (Ser313)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P41235
Dilution Range	Flow-Cyt=2ug/test,ICC=1:100
Expiration Date	12 months from date of receipt.