

## Product Datasheet

# AMPK alpha 1 (phospho-Thr172) antibody (PE) (orb504795)

<b>Description</b>	AMPK alpha 1 (phospho-Thr172) antibody (PE)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat, Sheep
<b>Conjugation</b>	PE
<b>Tested Applications</b>	IF
<b>Immunogen</b>	KLH conjugated Synthesised phosphopeptide derived from human AMPK alpha 2 around the phosphorylation site of Thr172
<b>Target</b>	phospho-AMPK alpha-2 (Thr172)
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">P54646</a>
<b>Expiration Date</b>	12 months from date of receipt.