

Product Datasheet

SHP2 (phospho-Tyr580) antibody (PE) (orb504857)

Description	Rabbit polyclonal antibody against SHP2 (phospho-Tyr580) conjugated to PE
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE
Tested Applications	ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PTPN11/SHP2 around the phosphorylation site of Tyr580 RV(p-Y)E
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	66kDa
Uniprot ID	Q06124
Dilution Range	ICC: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.