

Product Datasheet

MSK1 (phospho-Ser376) antibody (PE) (orb504915)

Description	Rabbit polyclonal antibody against MSK1 (phospho-Ser376) conjugated to PE
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	PE
Tested Applications	FC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human MSK1 around the phosphorylation site of Ser376 GY(p-S)FV
Target	RPS6KA5
Preservatives	Lyophilized or Liquid
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	90kDa
Uniprot ID	O75582
Dilution Range	IF: 1:50-200
Expiration Date	12 months from date of receipt.