
Product Datasheet

**HDAC4 (phospho-Ser246) antibody (PE)
(orb505017)**

Description

Rabbit polyclonal antibody against

to PE

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

PE

Tested Applications

IF

Immunogen

KLH conjugated Synthesised phosphopeptide derived from human HDAC4 around the phosphorylation site of Ser246

Target

HDAC4

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

119kDa

Uniprot ID**P56524****Expiration Date**

12 months from date of receipt.