

Product Datasheet

AMPK Beta 1 + AMPK Beta 2 (phospho-Ser182 antibody (PE) (orb505034)

Description	AMPK Beta 1 + AMPK Beta 2 (phospho-Ser182 antibody (PE))
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	PE
Tested Applications	FC, IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human AMPK beta 1 around the phosphorylation site of Ser182
Target	phospho-AMPK beta 1 + AMPK beta 2 (Ser182/Ser184)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y478
Dilution Range	Flow-Cyt=2ug/Test,IF=1:100-500
Expiration Date	12 months from date of receipt.