

Product Datasheet

CaMK2 alpha (phospho-Tyr231) antibody (PE) (orb505606)

| | |
|----------------------------|--|
| Description | CaMK2 alpha (phospho-Tyr231) antibody (PE) |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | PE |
| Tested Applications | FC, IF |
| Immunogen | KLH conjugated Synthesised phosphopeptide derived from human CaMK2 alpha around the phosphorylation site of Tyr231 |
| Target | phospho-CaMK2 alpha (Tyr230) |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q9UQM7 |
| Dilution Range | Flow-Cyt=2ug/Test,IF=1:100-500 |
| Expiration Date | 12 months from date of receipt. |