
Product Datasheet

**HER3 (phospho-Tyr1197) antibody (PE)
(orb505650)**

Description

Rabbit polyclonal antibody against

PE

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

PE

Tested Applications

IF, IHC-Fr, IHC-P

Immunogen

KLH conjugated Synthesised phosphopeptide derived from human HER3 around the phosphorylation site of Tyr1289 EE(p-Y)EY

Target

ERBB3

Preservatives

Lyophilized or Liquid

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

148kDa

CAS Number

3

Uniprot ID**P21860****Dilution Range**

IF: 1:50-200

Expiration Date

12 months from date of receipt.