

Product Datasheet

Lyn (phospho-Tyr507) antibody (PE) (orb505667)

Description	Rabbit polyclonal antibody against Lyn (phospho-Tyr507) conjugated to PE
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE
Tested Applications	ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Lyn around the phosphorylation site of Tyr507 GQ(p-Y)QQ
Target	LYN
Preservatives	Lyophilized or Liquid
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	57kDa
Uniprot ID	P07948
Dilution Range	IF: 1:50-200
Expiration Date	12 months from date of receipt.