

Product Datasheet

MEK1(phospho-Thr286) antibody (PE) (orb505677)

Description	Rabbit polyclonal antibody against MEK1(phospho-Thr286) conjugated to PE
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	PE
Tested Applications	IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human MEK1 around the phosphorylation site of Thr286
Target	MAP2K2
Preservatives	0.01M TBS(pH7.4) with 0.03% Proclin300.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at 2-8 °C for one year. Protect from light. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	43kDa
Uniprot ID	P36507
Expiration Date	12 months from date of receipt.