

Product Datasheet

TAK1 (phospho-Thr184/187) antibody (PE) (orb505738)

Description	TAK1 (phospho-Thr184/187) antibody (PE)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE
Tested Applications	IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human TAK1 around the phosphorylation site of Thr184/187
Target	Phospho-TAK1(Thr184 + Thr187)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O43318
Expiration Date	12 months from date of receipt.