
Product Datasheet

**CDC25C (phospho-Thr48) antibody (PE)
(orb505751)**

Description

Rabbit polyclonal antibody against

	PE
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE
Tested Applications	ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human cdc25C around the phosphorylation site of Thr48
Target	CDC25C
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	53kDa
Uniprot ID	P30307
Dilution Range	IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.