


## Product Datasheet

### IRAK1 (phospho-Thr209) antibody (orb506111)

<b>Description</b>	IRAK1 (phospho-Thr209) antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, IHC-P
<b>Immunogen</b>	KLH conjugated Synthesised phosphopeptide derived from human IRAK1 around the phosphorylation site of Thr209 (RG(p-T)HN)
<b>Target</b>	Phospho-IRAK1 (Thr209)
<b>Preservatives</b>	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	78 kDa
<b>Uniprot ID</b>	<a href="#">Q62406</a>
<b>Dilution Range</b>	IHC-P=1:100-500, ICC=1:100, IHC-F=1:100-500, IF=1:100-500
<b>Expiration Date</b>	12 months from date of receipt.



ICC analysis of HeLa cell using Phospho-...