

## Product Datasheet

# Human IL-17A/IL-4 Dual Fluorospot Set (orb506692)

<b>Description</b>	Biorbyt Dual Fluorospot Sets include capture antibody for cytokine 1, FITC-conjugated detection anti...
<b>Reactivity</b>	Human
<b>Note</b>	For research use only
<b>Application notes</b>	To receive kit plus non-sterile plates (MSIPN), please add 'P' at the end of the Cat N°. To receive kit plus sterile plates (MSIPS), please add 'S' at the end of the Cat N°.
<b>Expiration Date</b>	Please enquire.