

## Product Datasheet

### DIM-C-pPhOCH3(C-DIM5) (orb507407)

<b>Description</b>	C-DIM5 is a Nur77 agonist. Activation of the orphan nuclear receptor Nur77 by C-DIM5 is associated w...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{24}H_{20}N_2O$
<b>MW</b>	352.4284
<b>CAS Number</b>	[33985-68-1]
<b>SMILES</b>	<chem>COC1=CC=C(C(C2=CNC3=C2C=CC=C3)C4=CNC5=C4C=CC=C5)C=C1</chem>
<b>Expiration Date</b>	12 months from date of receipt.